



**A revolution in ultrasound-based diagnosis**

*The first test for early detection  
of preterm birth*



## AN UNSOLVED CLINICAL NEED

Prematurity \* is the leading cause of death for children and newborns in developed countries, the second worldwide, and presents an approximate global prevalence of 10%. Each year more than 15 million of preterm babies that can suffer complications due to this condition are born globally. Prematurity contributes to long-term growth impairment and significant long-term morbidity such as cognitive, visual and learning impairments. More than 50% of the total of preterm births are spontaneous preterm births (sPTB): unpredictable situations beyond a routine follow-up of pregnancy. Prematurity prevalence has remained stable during the last decades because the current standard test for its detection is based on the measurement of the cervical length (CL) \*\*. CL measurement performs very low in general population. For this reason, there is an open debate about its application between the guidelines which are defined among the different scientific societies (SMFM, ACOG and ISUOG).

### quantusPREMATURITY: the first fully automatic test for early detection of spontaneous preterm birth (sPTB)

- **Early detection:** quantusPREMATURITY is the first test that allows an early detection of the risk of sPTB before cervical shortening has already occurred through the automatic analysis of an ultrasound image of the cervix. quantusPREMATURITY brings the opportunity to avoid the variability when measuring cervical length, simplifying this process, and being more efficient thanks to the advanced image processing techniques and the Artificial Intelligence (AI).
- **Multi-Result:** quantusPREMATURITY allows sPTB prediction at two levels of pregnancy (before 34 and 37 weeks) for being able to foresee different scenarios of pregnancy's monitoring.
- **Reliable:** quantusPREMATURITY results presented better performance and greater strength than those test based on the CL.

#### Comparison of quantusPREMATURITY with the Cervical Length (CL):

|             | CL <sup>^</sup> | quantusPREMATURITY <sup>^</sup> | CL <sup>^</sup> | quantusPREMATURITY <sup>^</sup> |
|-------------|-----------------|---------------------------------|-----------------|---------------------------------|
| Sensitivity | 10.3%           | 16.5%                           | 10.0%           | 35.2%                           |
| Specificity | 99.1%           | 99.6%                           | 98.8%           | 99.7%                           |
| PPV         | 32.9%           | 52.2%                           | 10.0%           | 50.7%                           |
| NPV         | 95.7%           | 95.9%                           | 98.4%           | 98.9%                           |

<sup>^</sup> Average values in trials (references 9-10)

\* Defined as a condition suffered by babies born before 37 weeks of gestational age.

\*\* Transvaginal ultrasound-based technique which allows to evaluate cervix characteristics (uterine cervix) by measuring its length to determine the risk of preterm birth. Also known as cervical length.

## HOW TO USE quantusPREMATURITY?

Using quantusPREMATURITY is very easy, only 3 simple steps are required:




1. Acquire an ultrasound image

2. Upload the image to the web application

3. Receive the results report

### Step 1: Acquire an ultrasound image

quantusPREMATURITY requires a transvaginal ultrasound image of the uterine cervix in DICOM format. A simple guide that shows how to make acquisitions is available in the application.





**Step 2:** Use the web application of quantusPREMATURITY to analyze the image

This application is a simple tool that allows you to send the image that you want to analyze.

To do this, just follow these four simple steps and complete the analysis.



Upload

The DICOM image.  
More than one image  
according to convenience  
are allowed



Identify

The clinical information



Select

The desired image to  
be analyzed



Send

The sample  
to be analyzed

**Step 3:** Obtain the result from the application in few minutes

**Spontaneous Preterm Birth Test**

**quantusPREMATURITY**

**YOUR LOGO HERE**

**Antenatal & Provider Information**

Patient Name: [blank] Clinic Name: [blank]  
 Address: [blank] Center Name: [blank]  
 Patient ID: [blank] Referring Physician/Consultant: [blank]  
 Specimen: [blank] Doctor Name: [blank]  
 Date of Birth (MM/DD/YYYY): [blank] Referral Date: [blank]  
 Sample ID: [blank] Date Received: 12/11/2018

**Sample Information**

USACID NUMBER (REF): [blank]  
 Referring Center: [blank]  
 REPORT DATE: 12/11/2018 13:00  
 Referring Doctor: [blank]

**Test Result**

before 34 WEEKS' GESTATION:

RESULT: Low risk  
 Risk given by quantusPREMATURITY / Baseline Risk(%)  
 1.1 % / 1-2%

before 37 WEEKS' GESTATION:

RESULT: Low risk  
 Risk given by quantusPREMATURITY / Baseline Risk(%)  
 4.1 % / 4-10%

**SPONTANEOUS PRETERM BIRTH RISK**

**Additional reports**

**TRANSUR**  
 Technical Response  
 Manual Print Storage

## ¿WHEN TO USE quantusPREMATURITY?

During the 2nd trimester ultrasound-based morphological screening, the transvaginal scan to assess the risk of sPTE may be useful to estimate the need to carry out additional checks or launch specific therapeutic strategies for reducing this risk. sPTB risk prediction before the week 34 is especially relevant at clinical level, since it represents the gestational age with the largest number of fetal complications associated with prematurity.

quantusPREMATURITY is indicated for single pregnancies between 19 and 24 weeks of gestation to determine the risk of sPTB as a consequence of a premature cervical remodelling. It is designed to be used as a tool for universal screening in general population, regardless of presenting risk factors associated with preterm birth.

However, those patients with increased risk of sPTE WITH previous sPTB, background of conization and/or presence of uterine malformation, could particularly benefit from this test in order to adapt its follow-up or treatment in a personalized way.

For example, during ultrasound scanning of a first pregnancy of 20 weeks gestation, where measurement of cervical length is 25 mm and directed anamnesis, the patient does not present any other risk factor for preterm birth, the basal risk for sPTB < 34 weeks is 1.7% and < 37 weeks is 4.8%. A transvaginal ultrasound to evaluate the risk of PPE, together with the morphological ultrasound, is performed. A result of 'low risk' using quantusPREMATURITY would reduce the risk of sPTB < 34 weeks at 1.1% and < 37 weeks at 4.1%, while a 'high risk' result would increase the risk at 51.6% for sPTB < 34 weeks and a 73.3% for sPTB < 37 weeks. Therefore, early knowing the risk of sPTB by using a simple transvaginal ultrasound can clearly.





## quantusPREMATURITY OFFERS AN EXPERIENCE WITHOUT LIMITS

- ✓ Unrestricted 24/7 access: only having internet connection, you can use quantusPREMATURITY and check the results at any time and from anywhere.
- ✓ No installation is required: quantusPREMATURITY has been designed thinking about a simple initial use and does not require any software download or installation.
- ✓ Great compatibility: quantusPREMATURITY is compatible with most of the web browsers, as well as with the most common ultrasound machines used in obstetrics and Gynaecology.

## quantusPREMATURITY OFFERS A HIGH ECONOMIC VALUE

- ✓ Does not require any initial investment in infrastructure.
- ✓ Pay per use: pay only per analysis. Without mandatory minimum consumption or fee of monthly maintenance.
- ✓ 15 days freeTest available without commitment for adapting our technology to your clinical practice.

For further information please contact us at:  
[sales@transmuralbiotech.com](mailto:sales@transmuralbiotech.com)

## WHY DOES quantusPREMATURITY WORK?

Changes that occur at histological level in a tissue (proportion of collagen, fat or water, among others) affect the dispersion of ultrasonic signals. These signals constitute the basis for generation of ultrasound images. Quantitative analysis of ultrasound images, together with the advanced techniques based on Artificial Intelligence, allows to detect very subtle changes in tissues that are not visible to the human eye. These changes could be relevant information from tissue microstructure related to a pathology.

The uterine cervix is a potential candidate for the application of this type of technology, since a cervical remodelling must occur during pregnancy. Cervical remodelling process consists of microstructural changes of the cervix which starts within the first trimester of pregnancy and lasts until term. An early or premature cervical remodelling may precede sPTB. Therefore, an early identification of a premature cervical remodelling could reduce potential maternal and perinatal complications associated with sPTB.

quantusPREMATURITY provides an alternative to the early detection for the prediction of the sPTB risk reliably; reaching levels of precision, performance and reproducibility without precedent in the state of the art.

### References:

- [1] Chawanpaiboon S, Vogel JP, Moller AB, Lumbiganon P, Petzold M, Hogan D, Landoulsi S, Jampathong N, Kongwattanakul K, Laopaiboon M, Lewis C. Global, regional, and national estimates of levels of preterm birth in 2014: a systematic review and modelling analysis. *The Lancet Global Health*. 2018 Oct 30.
- [2] Romero R, Espinoza J, Kusanovic JP, Gotsch F, Hassan S, Erez O, Chaiworapongsa T, Mazor M. The preterm partition syndrome. *BJOG: An International Journal of Obstetrics&Gynaecology*. 2006 Dec 1;113(s3):17-42.
- [3] Mayers KM, Feltovich H, Mazza E, Vink J, Bajka M, Wapner RJ, Hall TJ, House M. The mechanical role of the cervix in pregnancy. *Journal of biomechanics*. 2015 Jun 25;48(9):1511-23.
- [4] Temming LA, Durst JK, Tuuli MG, Stout MJ, Dicke JM, Macones GA, Cahill AG. Universal cervical length screening: implementation and outcomes. *American Journal of Obstetrics&Gynecology*. 2016 Apr 1;214(4):523-e1.
- [5] Son M, Grobman WA, Ayala NK, Miller ES. A universal mid-trimester transvaginal cervical length screening program and its associated reduced preterm birth rate. *American Journal of Obstetrics&Gynecology*. 2016 Mar 1;214(3):365-e1.
- [6] Orzechowski KM, Boelig RC, Baxter JK, Berghella V. A universal transvaginal cervical length screening program for preterm birth prevention. *Obstetrics&Gynecology*. 2014 Sep 1;124(3):520-5.
- [7] Khalil MR, Thorsen P, Uldbjerg N. Cervical ultrasound elastography may hold potential to predict risk of preterm birth. *Dan Med J*. 2013 Jan 1;60(1):A4570.
- [8] Baños N, Murillo-Bravo C, Julià C, Migliorelli F, Perez-Moreno A, Rios J, Gratacos E, Valentin L, Palacio M. Mid-trimester sonographic cervical consistency index to predict spontaneous preterm birth in a low-risk population. *Ultrasound in Obstetrics&Gynecology*. 2018 May;51(5):629-36.
- [9] Baños N, Perez-Moreno A, Migliorelli F, Triginer L, Cobo T, Bonet-Carne E, Gratacos E, Palacio M. Quantitative analysis of the cervical texture by ultrasound and correlation with gestational age. *Fetal diagnosis and therapy*. 2017;41(4):265-72.
- [10] Baños N, Perez-Moreno A, Julià C, Murillo-Bravo C, Coronado D, Gratacos E, Deprest J, Palacio M. Quantitative analysis of cervical texture by ultrasound in mid-pregnancy and association with spontaneous preterm birth. *Ultrasound in Obstetrics&Gynecology*. 2018 May;51(5):637-43.



[www.quantusprematurity.com](http://www.quantusprematurity.com)



ROBUST



RELIABLE



FOR ALL WOMEN



Transmural Biotech S.L., CIF: B65084675.

C/ Sabino Arana, n38, 1º 1ª, 08028 Barcelona. Spain

Revisión 1 12.03.2018

quantusPREMATURETY is a medical device according  
to the EU regulations for CE mark