



***REVOLUTIONIZING THE
DETECTION OF MALIGNANT
DERMATOLOGICAL LESIONS***

- ✓ 100% non-invasive
- ✓ Reliable results in 10 minutes
- ✓ Excellent sensitivity and specificity values



AN UNRESOLVED CLINICAL NEED

- Malignant dermatologic lesions have increased significantly in incidence worldwide in white populations.
- If detected early, the lesion has a high chance of being treated.
- However, this has been exclusively performed on people with risk factors and there is no convincing evidence on the follow-up of low-risk groups of patients.

QUANTUSSKIN - ANALYSIS AND CLASSIFICATION OF DERMATOSCOPIC IMAGES FOR MALIGNANCY RISK ASSESSMENT

- quantusSKIN is a software that analyzes and classifies dermoscopic images to determinate the risk of malignancy of skin lesions.
- Non-invasive: quantusSKIN is a non-invasive test that predicts the risk of malignancy of different skin lesions through a photograph or a dermoscopic image.
- Fast: quantusSKIN generates accurate results in few minutes.

| Sensitivity | Specificity | PPV | NPV |
|-------------|-------------|-------|-------|
| 89,6% | 85,2% | 52,6% | 97,8% |

* Sensitivity: Proportion of negative cases correctly identified by the algorithm. It is the number of items correctly identified as negative out of the total number of negatives.

* Specificity: Proportion of positive cases correctly identified by the algorithm. It is the number of items correctly identified as positive over the actual total number of positives.

* PPV: Positive Predictive Value.

* NPV: Negative Predictive Value.


HOW TO USE quantusSKIN

Using quantusSKIN is simple, it only requires 3 steps:



Step1: ACQUIRE A DERMOSCOPIC IMAGE

quantusSKIN requires a skin image in JPG or PNG format captured through a smartphone, reflex camera or similar, always free of active light filters. A dermatoscope can also be utilized when magnification markers or size markers are not used. The application provides a simple guide that explains how to proceed with the acquisitions.



Step 2: Use the quantusSKIN medical app to analyze the image

The application allows the user to send the image that wants to analyze by following three simple steps.



Upload

Upload the
JPG or PNG
image



Select


Select the
desired image





Send

Send the
sample for
analysis

Step 3: Obtain the result of the application within a few minutes.

Screening test for malignant skin lesion 

| Patient & Provider Information | |
|--------------------------------|---------------------------|
| PATIENT NAME: | CLINIC NAME: |
| patient name | Transmural Biotech |
| PATIENT ID: | REFERRING/ORDER CLINICAL: |
| 122345 | David Fernández |
| quantusSKIN ID: | REPORT DATE: |
| btech-1084 | (dd/mm/yyyy) 11/01/2021 |

| Sample Information | Test Result quantusSKIN |
|---|--|
|  | quantusSKIN ID: |
| VS ACQUISITION DATE: (dd/mm/yyyy) 04/01/2021 | btech-1084 |
| REQUEST DATE: (dd/mm/yyyy hh:mm) 11/01/2021 15:58 | quantusSKIN Risk: |
| | 10% 0% 25% 50% 75% 100% |
| | Authorized by  |

TEST DESCRIPTION:
quantusSKIN® offers an automatic assessment of the skin cancer risk using the quantitative texture analysis of a dermoscopic image of the skin lesion. Quality of image and acquisition is relevant and must be taken following quantusSKIN requirements, quantusSKIN has been developed for clinical research only.

quantusSKIN is considered for clinical research use only and not for clinical use. quantusSKIN has been developed by Transmural Biotech, Biotecnología, S.L. (Barcelona, ES, since 1998). Biotechnología Spain is registered with the Spanish regulatory authorities for the purpose of clinical studies and reports to the corresponding Spanish authorities. Transmural Biotech is the in charge of processing personal data in order to offer you medical treatment. You can exercise the rights of access, rectification, cancellation, opposition, request, data portability, report and restriction of processing consulting us at info@transmuralbiotech.com

WHEN TO USE quantusSKIN

quantusSKIN is a non-invasive, fast and easy-to-use test that allows the detection of malignant dermatological lesions through dermoscopic images. Its technology is based on quantitative analysis of dermoscopic image texture. By simply analyzing and classifying images, quantusSKIN determines in few minutes the probability of a skin lesion.

quantusSKIN design has been focused on general population with the purpose of being a tool for the detection of malignant skin lesions such as melanoma, basal cell carcinoma or squamous cell carcinoma. Moreover, it allows the screening of patients with risk factors and the prioritization in waiting lists.

quantusSKIN classifies skin lesions in benign or malignant without the need of or in addition to visual inspection from a specialist via a dermatoscope. The specialist, classifies the images by visual patterns and quantusSKIN gives a percentage risk of malignancy.



AN INNOVATIVE MEDICAL SOLUTION:

- ✓ **Unrestricted 24-hour access:** It is essential to have an internet connection to use quantusSKIN and review results at any time and from any location.
- ✓ **No installation required:** quantusSKIN has been designed in such a way that its initial use is simple since it does not require the download or installation of any software.
- ✓ **High compatibility:** quantusSKIN is compatible with most browsers. The model can be used for web-based as well as primary devices.

quantusSKIN OFFERS HIGH ECONOMIC VALUE:

- ✓ **NO initial investment in infrastructure required!**
- ✓ **Pay-as-you-go:** Pay only for each test you order!
- ✓ **FREE 30-day trial available, no obligation!**
- ✓ **Add more value to your clinic and increase your profits!**



To get a FREE 30-day trial,
please contact us at
sales@transmuralbiotech.com



+34 931 190 929



+34 626 667 989

WHY DOES quantusSKIN WORK?

An automated support tool requires minimal or no physician intervention to obtain a result. Over the past few years, research has been focused on automated algorithms to improve current imaging-based clinical diagnosis. The rise of Artificial Intelligence techniques, and especially Deep Learning, has increased the number of studies using this type of algorithm in diagnostic dermatology.

Recently published studies show that skin cancer detection using trained Deep Learning models can achieve high accuracy in diverse populations.

It provides quantitative comparisons on how model performance can vary across datasets consisting of malignant dermatological lesions of different disease severity and ethnicity.

quantusSKIN is presented as a novel Artificial Intelligence method based on state-of-the-art Deep Learning. Several studies prove the efficacy of the analysis method in the quantitative analysis method proposed by quantusSKIN.

The technology is based on performing a quantitative analysis of the texture of the cutaneous Nevus image obtained using a smartphone, reflex camera or dermatoscope. This analysis allows to identify patterns associated with specific pathologies and to determine the risk of malignancy of the skin lesion. According to the literature, the various tests and tools used by the dermatologist give an individual sensitivity of 75-84% (Dermatol,2008)⁹. While quantusSKIN has obtained in its tests a sensitivity of 85.6% (Coronado-Gutiérrez, et al.,2021)¹⁶.

References

- [1] U. Leiter, T. Eigentler, and C. Garbe, "Epidemiology of Skin Cancer BT - Sunlight, Vitamin D and Skin Cancer," in *Advances in experimental medicine and biology*, vol. 810, J. Reichrath, Ed. Springer New York, 2014, pp. 120–140.
- [2] C. Garbe and U. Leiter, "Melanoma epidemiology and trends," *Clin. Dermatol.*, vol. 27, no. 1, pp. 3–9, Jan. 2009, doi: 10.1016/j.clindermatol.2008.09.001.
- [3] G. P. Guy, S. R. Machlin, D. U. Ekwueme, and K. R. Yabroff, "Prevalence and costs of skin cancer treatment in the U.S., 2002–2006 and 2007–2011," *Am. J. Prev. Med.*, vol. 48, no. 2, pp. 183–187, Feb. 2015, doi: 10.1016/j.amepre.2014.08.036.
- [4] H. W. Rogers, M. A. Weinstock, S. R. Feldman, and B. M. Coldiron, "Incidence estimate of nonmelanoma skin cancer (keratinocyte carcinomas) in the us population, 2012," *JAMA Dermatology*, vol. 151, no. 10, pp. 1081–1086, Oct. 2015, doi: 10.1001/jamadermatol.2015.1187.
- [5] R. L. Siegel, K. D. Miller, and A. Jemal, "Cancer statistics, 2017," *CA. Cancer J. Clin.*, vol. 67, no. 1, pp. 7–30, Jan. 2017, doi: 10.3322/caac.21387.
- [6] "Melanoma Warning Signs and Images - The Skin Cancer Foundation." <https://www.skincancer.org/skin-cancer-information/melanoma/melanoma-warning-signs-and-images/> (accessed Sep. 29, 2020).
- [7] H. A. Haenssle et al., "Association of patient risk factors and frequency of nevus-associated cutaneous melanomas," *JAMA Dermatology*, vol. 152, no. 3, pp. 291–298, Mar. 2016, doi: 10.1001/jamadermatol.2015.3775.
- [8] P. Tschandi and P. Doz Philipp Tschandi, "Sequential digital dermatoscopic imaging of patients with multiple atypical nevi," *Rev. J Dermatol Pr. Concept*, vol. 8, no. 3, pp. 231–237, 2018, doi: 10.5826/dpc.0803a16.
- [9] M. E. Vestergaard, P. Macaskill, P. E. Holt, and S. W. Menzies, "Dermoscopy compared with naked eye examination for the diagnosis of primary melanoma: A meta-analysis of studies performed in a clinical setting," *Br. J. Dermatol.*, vol. 159, no. 3, pp. 669–676, Sep. 2008, doi: 10.1111/j.1365-2133.2008.08713.x.
- [10] H. Kittler, H. Pehamberger, K. Wolff, and M. Binder, "Diagnostic accuracy of dermoscopy," *Lancet Oncology*, vol. 3, no. 3, Lancet Publishing Group, pp. 159–165, Mar. 01, 2002, doi: 10.1016/S1470-2045(02)00679-4.
- [11] A. C. Geller, S. M. Swetter, K. Brooks, M. F. Demierre, and A. L. Yaroch, "Screening, early detection, and trends for melanoma: Current status (2000–2006) and future directions," *Journal of the American Academy of Dermatology*, vol. 57, no. 4, Mosby, pp. 555–572, Oct. 01, 2007, doi: 10.1016/j.jaad.2007.06.032.
- [12] A. Rosenberg and J. H. Meyerle, "Total-body photography in skin cancer screening: The clinical utility of standardized imaging," *Cutis*, vol. 99, no. 5, pp. 312–316, May 2017, Accessed: Sep. 29, 2020. [Online]. Available: <https://europepmc.org/article/med/28632800>.
- [13] N. C. F. Codella et al., "Skin lesion analysis toward melanoma detection: A challenge at the 2017 International symposium on biomedical imaging (ISBI), hosted by the international skin imaging collaboration (ISIC)," in *Proceedings - International Symposium on Biomedical Imaging*, May 2018, vol. 2018-April, pp. 168–172, doi: 10.1109/ISBI.2018.8363547.
- [14] P. Tschandi, C. Rosendahl, and H. Kittler, "Data descriptor: The HAM10000 dataset, a large collection of multi-source dermatoscopic images of common pigmented skin lesions," *Sci. Data*, vol. 5, no. 1, pp. 1–9, Aug. 2018, doi: 10.1038/sdata.2018.161.
- [15] M. Combalia et al., "BCN20000: Dermoscopic Lesions in the Wild," Aug. 2019, Accessed: Jul. 01, 2020. [Online]. Available: <http://arxiv.org/abs/1908>.
- [16] Coronado-Gutiérrez, D., López, C., & Burgos-Artizzu, X. (2021). Skin cancer high-risk patient screening from dermoscopic images via Artificial Intelligence: an online study. doi: 10.1101/2021.02.04.21251132



www.quantusSKIN.org



NON INVASIVE



RELIABLE



FAST



To get a **FREE 30-day trial**,
please contact us at
sales@transmuralbiotech.com



Tel: +34 931 190 929 ☎ +34 626 667 989

Transmural Biotech S.L, CIF: B65084675.

S/ Beethoven 15, Floor 4 Office 18 08021 Barcelona, Spain

Email.: sales@transmuralbiotech.com

📷 quantustb

🌐 TRANSMURAL BIOTECH



Revision 4 - 06 . 07 . 2021