



REVOLUTIONIZING DIAGNOSTIC IMAGING

The first test for the early detection of spontaneous preterm birth

AN UNMET CLINICAL NEED



Prematurity is the leading cause of infant mortality in developed countries, the second around the world and has a global prevalence of approximately 10%. Every year more than 15 million premature babies are born in the world who can suffer complications due to this condition. Prematurity contributes in the long term to suffering from growth disability and other pathologies such as cognitive morbidity among others.

The incidence of prematurity has remained stable over the past few decades, as the current reference test for detection is based on the measurement of cervical length (CL).

quantusPREMATURITY ITY The first test

- **Early detection:** quantusPREMATURITY is the first test that allows an early detection of the risk of sPTB before the cervix is shortened, through the automatic analysis of an ultrasound of the uterine cervix. It provides the opportunity to avoid cervical length variability by simplifying the process, making it more efficient thanks to advanced image processing techniques and Artificial Intelligence (AI).
- **Multi-result:** With quantusPREMATURITY we obtain the prediction of sPTB at two levels of gestation (before 34 and 47 weeks) in order to predict different follow-up scenarios in pregnancy.
- **Reliable:** quantusPREMATURITY results are better performing and more robust than CL-based results.

Comparison of quantusPREMATURITY with Cervicometry (CL):

Risk	Test	Sensitivity	Specificity	PPV	NPV
PPE < 34 Weeks		32,4%	99,5%	53,7%	98,8%
	CL	10,0%	98,8%	10,0%	98,4%
PPE < 37 Weeks		13,2%	98,7%	36,8%	95,8%
	CL	10,3%	99,1%	32,9%	95,7%

*Sensitivity: probability that the test will identify as sick the one who is actually sick.

*Specificity: probability that the test will identify as not sick the one who is not actually sick.

*PPV and NPV: Positive Predictive Value and Negative Predictive Value.

HOW TO USE QUANTUSPREMATURITY?

Using quantusPREMATURITY is easy, it only requires 3 simple steps:



**1. Purchase
one
echography**




**2. Upload the
image
to the web
application**



**3. Get the
results**

Paso 1: Getting an ultrasound quantusPREMATURITY requires a transvaginal ultrasound of the neck in DICOM format. A simple guide is available in the application showing how to make acquisitions



Step 2:

Use the **quantusPREMATURITY** web application to analyze the image.

This application is a simple tool that allows you to send to the system the image you want to analyze. To do this, you only have to follow three simple steps to complete the analysis:



Subir

The DICOM image
Allowed to upload
more
of an image
according to your
convenience



Seleccionar

The desired image
to be analyzed



Enviar

The sample to be
analyzed

Step 3: Get the result of the application in few minutes.

Spontaneous Preterm Birth Test



www.quantuspremat.com

Patient & Provider Information	
PATIENT NAME:	CLINIC NAME:
Patient Name	Transmural Biotech
PATIENT ID:	ALTERNATIVE CLINIC:
PatientID	Dr. Name
quantusPREMATURITY ID:	REPORT DATE:
TRANS-06	MM/DD/YYYY 10/01/2019

Sample Information	Test Result
	before 34 WEEKS' GESTATION: RESULT: Low risk Risk given by quantusPREM(*) / Baseline Risk(**): 1.1 % / 1-2%
ON ACQUISITION DATE: MM/DD/YYYY	before 37 WEEKS' GESTATION: RESULT: High risk Risk given by quantusPREM(*) / Baseline Risk(**): 46.7 % / 4-10%
REQUEST DATE: MM/DD/YYYY (H:MM)	

SPONTANEOUS PRETERM BIRTH RISK	Authorized signatory
	 Technical Responsible: Álvaro Pérez Moreno

(*) Comparison may given by quantusPREM(**) based on special population of low-risk and high-risk women with singleton pregnancies.
(**) Baseline risk has been established from previous large studies of spontaneous Preterm Birth (SPB) in general populations. Please be aware that baseline SPB risk is high. Based on an overall SPB risk of 10%, gestation 34 and other cases, risk of SPB is 10%.

quantusPREMATURITY offers an automatic assessment of the spontaneous Preterm Birth (SPB) risk using the quantitative analysis of a cervical image. Quality of image and acquisition is relevant and must be taken. Primary contact: quantuspremat@transmural.com. It has been validated in a general population of singleton pregnancies, women of singleton pregnancies from 32 to 34 weeks. quantusPREM(PT) is not intended as a diagnostic tool for identifying an additional information to the clinician in the management of the patient.

quantusPREMATURITY is a service provided by Transmural Biotech S.L. (Spain) and is not intended for use in any other country. The use of this service is subject to the terms and conditions of the service provided by Transmural Biotech S.L. For more information, please contact: info@transmural.com or +34 91 540 00 00. The use of this service is subject to the terms and conditions of the service provided by Transmural Biotech S.L. For more information, please contact: info@transmural.com or +34 91 540 00 00.

WHEN TO USE QUANTUSPREMATURIT

During morphological screening ultrasound in the second trimester, it may be very useful to associate transvaginal examination to assess the risk of spontaneous preterm birth and to estimate the need for additional controls or initiate specific therapeutic strategies, in order to reduce this risk. The prediction of the risk of EPP before week 34 is especially relevant at the clinical level, since it represents the gestational age with the highest number of fetal complications associated with prematurity.

quantusPREMATURITY is indicated for single pregnancies between 19 and 24 weeks gestation to determine the risk of spontaneous preterm birth as a result of premature cervical remodeling. For those patients with an increased risk of EPP with a history of PEP in previous gestation, a history of conization and/or presence of uterine malformation, they could especially benefit from this test in order to adapt their follow-up and / or treatment in a personalized way.

For example, in a 20-week primigestation of 20 weeks, where the measurement of cervical length is 25 mm and in the directed anamnesis, the patient does not represent any other risk factor for preterm birth, the baseline risk for EPP for <34 weeks is 1.7% and for <37 weeks 4.8%. The patient undergoes, together with morphological ultrasound, a transvaginal ultrasound to assess the risk of EPP. With a "low risk" result in quantusPREMATURITY would reduce the risk of PPE in <34 weeks to 1.1% and in <37 weeks by 4.1% while a "high risk" result would increase it to 25.8% for <34 weeks and 46.7% for <37 weeks. Therefore, knowing the risk of EPP early through a simple transvaginal ultrasound can have a clear impact on the clinical management of the case to assess subsequent controls or start of personalized treatments.



AN EXPERIENCE WITHOUT LIMITS

- ✓ **24-hour unrestricted access:** Through an internet connection it is possible to use **quantusPREMATURITY** and review the results at any time and from anywhere.
No installation required: or download of any type of Software.
- ✓ **Great compatibility:** **quantusPREMATURITY** is compatible with most web browsers as well as devices used in medical practice.
- ✓

quantusFLM OFFERS A HIGH ECONOMIC VALUE

- ✓ **It does NOT require initial investment in infrastructure!**
Pay as you go: Pay only for each analysis you request!
- ✓
- ✓ **It brings more value to your clinic and increases your profits!**



WHY DOES QUANTUSPREMATURITY work?


The changes that occur at the histological level in a tissue, including the proportion of collagen, fat or water, among others, affect the dispersion of ultrasound signals. These signals form the basis for ultrasound imaging.

The quantitative analysis of ultrasound images, together with advanced Artificial Intelligence techniques, allows to detect very subtle changes in tissues that are not visible to the human eye. These changes could be relevant information of the microstructure of the tissue related to a pathology.

The uterine cervix is a potential candidate for the application of this type of technology, since it must be remodeled during pregnancy until the moment of delivery. The process of cervical remodeling consists of microstructural changes of the cervix that begin in the first trimester of pregnancy and end in childbirth.

Early or premature remodeling of the cervix may be related to PTB. Thus, early identification of a premature cervical remodeling process could reduce the possible maternal and perinatal complications associated with PTB.

quantusPREMATURITY provides an early detection alternative for reliable risk prediction of PTB; reaching levels of accuracy, performance and reproducibility unprecedented in the state of the art.






Precise



Effective



Fast

**If you want to try it
Contact us!**



Email.: sales@transmuralbiotech.com

 **quantustb**  **TRANSMURAL BIOTECH**

CE



**TRANSMURAL
BIOTECH**



+34 931 190



+34 626 667

Transmural Biotech S.L., CIF: B65084675.

C/ Beethoven 15 Planta 4 Desp. 18 08021 Barcelona